

## Product datasheet for **TB402858**

### p60 CAF1 (CHAF1B) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CHAF1B (NM_005441) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202858
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	p60 CAF1 (CHAF1B) Mouse Monoclonal Antibody [Clone ID: OTI4F7] (TA506840)
ACCN:	<a href="#">NM_005441</a> , <a href="#">NP_005432</a>
Synonyms:	CAF-1; CAF-IP60; CAF1; CAF1A; CAF1P60; MPHOSPH7; MPP7
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_005432</a>
Locus ID:	8208
Cytogenetics:	21q22.12-q22.13
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)