

Product datasheet for **TB402933**

CKS1 (CKS1B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CKS1B (NM_001826), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202933
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001826 , NP_001817
Synonyms:	CKS1; ckshs1; PNAS-16; PNAS-18
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001817
Locus ID:	1163
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Pathways in cancer, Small cell lung cancer



[View online »](#)