

Product datasheet for **TB402971**

ZKSCAN3 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZKSCAN3 (NM_024493) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202971
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_024493 , NP_077819
Synonyms:	dj874C20.1; dj874C20.1.; ZF47; zfp-47; Zfp47; ZFP306; ZNF306; ZNF309; ZSCAN13; ZSCAN35
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_077819
Locus ID:	80317
Cytogenetics:	6p22.1
Protein Families:	Transcription Factors



[View online »](#)