

## Product datasheet for **TB403148**

### BHMT CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BHMT (NM_001713) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC203148
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	BHMT Mouse Monoclonal Antibody [Clone ID: OTI3E11] (TA500959)
ACCN:	<u>NM_001713, NP_001704</u>
Synonyms:	BHMT1; HEL-S-61p
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001704</u>
Locus ID:	635
Cytogenetics:	5q14.1
Protein Pathways:	Cysteine and methionine metabolism, Glycine, serine and threonine metabolism, Metabolic pathways



[View online »](#)