

## Product datasheet for **TB403507**

### Seladin 1 (DHCR24) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DHCR24 (NM_014762) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC203507
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_014762</a> , <a href="#">NP_055577</a>
Synonyms:	DCE; Nbla03646; seladin-1; SELADIN1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_055577</a>
Locus ID:	1718
Cytogenetics:	1p32.3
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transmembrane
Protein Pathways:	Metabolic pathways, Steroid biosynthesis



[View online »](#)