

Product datasheet for **TB404051**

UQCRFS1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UQCRFS1 (NM_006003) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204051
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	UQCRFS1 Mouse Monoclonal Antibody [Clone ID: OTI4H8] (TA811461)
ACCN:	NM_006003 , NP_005994
Synonyms:	MC3DN10; RIP1; RIS1; RISP; UQCR5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005994
Locus ID:	7386
Cytogenetics:	19q12
Protein Pathways:	Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease



[View online »](#)