

Product datasheet for **TB404090**

EFEMP1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EFEMP1 (NM_001039349), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204090
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	EFEMP1 Mouse Monoclonal Antibody [Clone ID: OT11D9] (TA503772)
ACCN:	NM_001039349 , NP_001034438
Synonyms:	DHRD; DRAD; FBLN3; FBNL; FIBL-3; MLVT; MTLV; S1-5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001034438
Locus ID:	2202
Cytogenetics:	2p16.1
Protein Families:	Secreted Protein



[View online »](#)