

Product datasheet for **TB404349**

PMP22 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PMP22 (NM_153321), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204349
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PMP22 Mouse Monoclonal Antibody [Clone ID: OT15D1] (TA808964)
ACCN:	NM_153321 , NP_696996
Synonyms:	CIDP; CMT1A; CMT1E; DSS; GAS-3; GAS3; HMSNIA; HNPP; Sp110
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_696996
Locus ID:	5376
Cytogenetics:	17p12
Protein Families:	Transmembrane



[View online »](#)