

## Product datasheet for **TB404470**

### Cytokeratin 16 (KRT16) CytoSection

#### Product data:

|                                       |                                                                                                                                                                                                                         |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | CytoSections                                                                                                                                                                                                            |
| Description:                          | Transient overexpression of KRT16 (NM_005557) in HEK293T cells paraffin embedded controls for ICC/IHC staining                                                                                                          |
| Species:                              | Human                                                                                                                                                                                                                   |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                 |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC204470                                                                                                                                                                                                  |
| Tag:                                  | C-MYC/DDK                                                                                                                                                                                                               |
| Detection Antibodies:                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)                                                                                                                                                         |
| Target Detection Antibodies:          | Cytokeratin 16 (KRT16) Mouse Monoclonal Antibody [Clone ID: OTI7A4] (TA808338)                                                                                                                                          |
| ACCN:                                 | <a href="#">NM_005557</a> , <a href="#">NP_005548</a>                                                                                                                                                                   |
| Synonyms:                             | CK16; FNEPPK; K1CP; K16; KRT16A; NEPPK; PC1                                                                                                                                                                             |
| Storage:                              | Room Temperature, or 2-8°C for long term storage                                                                                                                                                                        |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                 |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                   |
| RefSeq:                               | <a href="#">NP_005548</a>                                                                                                                                                                                               |
| Locus ID:                             | 3868                                                                                                                                                                                                                    |
| Cytogenetics:                         | 17q21.2                                                                                                                                                                                                                 |



[View online »](#)