

Product datasheet for **TB404530**

ARHGEF1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ARHGEF1 (NM_004706), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204530
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_004706 , NP_004697
Synonyms:	GEF1; IMD62; LBCL2; LSC; P115-RHOGEF; SUB1.5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004697
Locus ID:	9138
Cytogenetics:	19q13.2
Protein Pathways:	Regulation of actin cytoskeleton, Vascular smooth muscle contraction



[View online »](#)