

Product datasheet for **TB404639**

FARSLB (FARSB) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FARSB (NM_005687) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204639
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	FARSLB (FARSB) Mouse Monoclonal Antibody [Clone ID: OTI4B3] (TA806739)
ACCN:	NM_005687 , NP_005678
Synonyms:	FARSLB; FRSB; HSPC173; NEDBLLA; PheHB; PheRS; RILDBC; RILDBC1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005678
Locus ID:	10056
Cytogenetics:	2q36.1
Protein Pathways:	Aminoacyl-tRNA biosynthesis



[View online »](#)