

Product datasheet for **TB404779**

CHMP5 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CHMP5 (NM_016410), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204779
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CHMP5 Mouse Monoclonal Antibody [Clone ID: OTI5C7] (TA808523)
ACCN:	NM_016410 , NP_057494
Synonyms:	C9orf83; CGI-34; HSPC177; PNAS-2; SNF7DC2; Vps60
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_057494
Locus ID:	51510
Cytogenetics:	9p13.3
Protein Pathways:	Endocytosis



[View online »](#)