

Product datasheet for **TB404861**

MTF1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MTF1 (NM_005955) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204861
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	MTF1 Mouse Monoclonal Antibody [Clone ID: OTI2F3] (TA805480)
ACCN:	NM_005955 , NP_005946
Synonyms:	MTF-1; ZRF
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005946
Locus ID:	4520
Cytogenetics:	1p34.3
Protein Families:	Transcription Factors



[View online »](#)