

Product datasheet for **TB404918**

LGALS3BP CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LGALS3BP (NM_005567) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204918
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	LGALS3BP Mouse Monoclonal Antibody [Clone ID: OTI3G8] (TA503455)
ACCN:	NM_005567 , NP_005558
Synonyms:	90K; BTBD17B; CyCAP; gp90; M2BP; MAC-2-BP; TANGO10B
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005558
Locus ID:	3959
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)