

Product datasheet for **TB404952**

XRCC1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of XRCC1 (NM_006297) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204952
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	XRCC1 Mouse Monoclonal Antibody [Clone ID: OTI2F8] (TA500880)
ACCN:	NM_006297 , NP_006288
Synonyms:	RCC; SCAR26
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006288
Locus ID:	7515
Cytogenetics:	19q13.31
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair



[View online »](#)