

Product datasheet for **TB405030**

UPRT CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UPRT (NM_145052), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205030
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UPRT Mouse Monoclonal Antibody [Clone ID: OT11A3] (TA506871)
ACCN:	NM_145052 , NP_659489
Synonyms:	FUR1; UPP
Storage:	RT
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_659489
Locus ID:	139596
Cytogenetics:	Xq13.3
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Pyrimidine metabolism



[View online »](#)