

## Product datasheet for **TB405066**

### Viperin (RSAD2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RSAD2 (NM_080657) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205066
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Viperin (RSAD2) Mouse Monoclonal Antibody [Clone ID: OTI4D12] (TA505799)
ACCN:	<u><a href="#">NM_080657</a></u> , <u><a href="#">NP_542388</a></u>
Synonyms:	cig5; cig33; vig1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_542388</a></u>
Locus ID:	91543
Cytogenetics:	2p25.2



[View online »](#)