

Product datasheet for **TB405070**

Aph 1b (APH1B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of APH1B (NM_031301), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205070
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_031301</u> , <u>NP_112591</u>
Synonyms:	APH-1B; PRO1328; PSFL; TAAV688
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_112591</u>
Locus ID:	83464
Cytogenetics:	15q22.2
Protein Families:	Transmembrane



[View online »](#)