

Product datasheet for **TB405137**

ERLIN1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ERLIN1 (NM_006459) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205137
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	ERLIN1 Mouse Monoclonal Antibody [Clone ID: OTI4B1] (TA809148)
ACCN:	NM_006459 , NP_006450
Synonyms:	C10orf69; Erlin-1; KE04; KEO4; SPFH1; SPG62
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006450
Locus ID:	10613
Cytogenetics:	10q24.31
Protein Families:	Transmembrane



[View online »](#)