

Product datasheet for **TB405236**

DPP8 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DPP8 (NM_197961), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205236
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DPP8 Mouse Monoclonal Antibody [Clone ID: OTI1D2] (TA504366)
ACCN:	NM_197961 , NP_932065
Synonyms:	DP8; DPRP-1; DPRP1; MST097; MSTP097; MSTP135; MSTP141
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_932065
Locus ID:	54878
Cytogenetics:	15q22.31
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)