

Product datasheet for **TB405289**

G protein alpha inhibitor 1 (GNAI1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GNAI1 (NM_002069) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205289
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	G protein alpha inhibitor 1 (GNAI1) Mouse Monoclonal Antibody [Clone ID: OTI2D1] (TA811907)
ACCN:	<u>NM_002069</u> , <u>NP_002060</u>
Synonyms:	Gi
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_002060</u>
Locus ID:	2770
Cytogenetics:	7q21.11
Protein Families:	Druggable Genome
Protein Pathways:	Axon guidance, Chemokine signaling pathway, Gap junction, Leukocyte transendothelial migration, Long-term depression, Melanogenesis, Progesterone-mediated oocyte maturation, Tight junction



[View online »](#)