

Product datasheet for **TB405441**

DPPA2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DPPA2 (NM_138815) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205441
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	DPPA2 Mouse Monoclonal Antibody [Clone ID: OTI1G10] (TA800361)
ACCN:	NM_138815 , NP_620170
Synonyms:	CT100; ECAT15-2; PESCRG1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_620170
Locus ID:	151871
Cytogenetics:	3q13.13
Protein Families:	Embryonic stem cells, Induced pluripotent stem cells, Stem cell - Pluripotency



[View online »](#)