

## Product datasheet for **TB405442**

### VTA1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of VTA1 (NM_016485) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205442
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_016485</a> , <a href="#">NP_057569</a>
Synonyms:	C6orf55; DRG-1; DRG1; HSPC228; LIP5; My012; SBP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_057569</a>
Locus ID:	51534
Cytogenetics:	6q24.1-q24.2
Protein Pathways:	Endocytosis



[View online »](#)