

Product datasheet for **TB405695**

Asialoglycoprotein Receptor 2 (ASGR2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ASGR2 (NM_080913), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205695
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Asialoglycoprotein Receptor 2 (ASGR2) Mouse Monoclonal Antibody [Clone ID: OTI2A12] (TA505099)
ACCN:	NM_080913 , NP_550435
Synonyms:	ASGP-R2; ASGPR2; CLEC4H2; HBXBP; HL-2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_550435
Locus ID:	433
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)