

Product datasheet for **TB405867**

SNAP47 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SNAP47 (NM_053052) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205867
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_053052 , NP_444280
Synonyms:	C1orf142; ESFI5812; HEL-S-290; HEL170; SNAP-47; SVAP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_444280
Locus ID:	116841
Cytogenetics:	1q42.13
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)