

Product datasheet for **TB405881**

beta II Tubulin (TUBB2A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TUBB2A (NM_001069) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205881
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	beta II Tubulin (TUBB2A) Mouse Monoclonal Antibody [Clone ID: OTI3E6] (TA506590)
ACCN:	NM_001069 , NP_001060
Synonyms:	CDCBM5; TUBB; TUBB2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001060
Locus ID:	7280
Cytogenetics:	6p25.2
Protein Families:	Druggable Genome
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection



[View online »](#)