

## Product datasheet for **TB405911**

### Glypican 3 (GPC3) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of GPC3 (NM_004484), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining   |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC205911  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| ACCN:                                 | <a href="#">NM_004484</a> , <a href="#">NP_004475</a>   |
| Synonyms:                             | DGSX; GTR2-2; MXR7; OCI-5; SDYS; SGB; SGBS; SGBS1   |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <a href="#">NP_004475</a>   |
| Locus ID:                             | 2719  |
| Cytogenetics:                         | Xq26.2  |
| Protein Families:                     | Druggable Genome  |



[View online »](#)