

Product datasheet for **TB406948**

TRAF5 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TRAF5 (NM_001033910), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206948
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TRAF5 Mouse Monoclonal Antibody [Clone ID: OTI2G3] (TA810290)
ACCN:	NM_001033910 , NP_001029082
Synonyms:	MGC:39780; RNF84
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001029082
Locus ID:	7188
Cytogenetics:	1q32.3
Protein Families:	Druggable Genome
Protein Pathways:	Pathways in cancer, Small cell lung cancer



[View online »](#)