

Product datasheet for **TB407051**

Kir2.2 (KCNJ12) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNJ12 (NM_021012) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207051
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_021012</u> , <u>NP_066292</u>
Synonyms:	hIRK; hIRK1; hkir2.2x; IRK-2; IRK2; kcnj12x; KCNJN1; Kir2.2; Kir2.2v
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_066292</u>
Locus ID:	3768
Cytogenetics:	17p11.2
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane



[View online »](#)