

## Product datasheet for **TB407345**

### **P2RY5 (LPAR6) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of LPAR6 (NM_005767), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC207345
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>ACCN:</b>	<u><a href="#">NM_005767</a></u> , <u><a href="#">NP_005758</a></u>
<b>Synonyms:</b>	ARWH1; HYPT8; LAH3; LPA-6; P2RY5; P2Y5
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u><a href="#">NP_005758</a></u>
<b>Locus ID:</b>	10161
<b>Cytogenetics:</b>	13q14.2
<b>Protein Families:</b>	Druggable Genome, GPCR, Transmembrane
<b>Protein Pathways:</b>	Neuroactive ligand-receptor interaction



[View online »](#)