

Product datasheet for **TB407360**

DARS2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DARS2 (NM_018122) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207360
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_018122 , NP_060592
Synonyms:	ASPRS; LBSL; MT-ASPRS; mtAsPRS
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060592
Locus ID:	55157
Cytogenetics:	1q25.1
Protein Pathways:	Aminoacyl-tRNA biosynthesis



[View online »](#)