

## Product datasheet for **TB407587**

### beta 2 Microglobulin (B2M) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of B2M (NM_004048) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207587
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OT14C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	beta 2 Microglobulin (B2M) Mouse Monoclonal Antibody [Clone ID: OT11A9] (TA802776)
ACCN:	<a href="#">NM_004048</a> , <a href="#">NP_004039</a>
Synonyms:	IMD43
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_004039</a>
Locus ID:	567
Cytogenetics:	15q21.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Antigen processing and presentation



[View online »](#)