

## Product datasheet for **TB407709**

### G CSF (CSF3) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CSF3 (NM_172220), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207709
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	G CSF (CSF3) Mouse Monoclonal Antibody [Clone ID: OT13F12] (TA810788)
ACCN:	<a href="#">NM_172220</a> , <a href="#">NP_757374</a>
Synonyms:	C17orf33; CSF3OS; GCSF
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_757374</a>
Locus ID:	1440
Cytogenetics:	17q21.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway



[View online »](#)