

Product datasheet for **TB407795**

HSPA6 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HSPA6 (NM_002155) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207795
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	HSPA6 Mouse Monoclonal Antibody [Clone ID: OTI1D12] (TA501950)
ACCN:	<u>NM_002155</u> , <u>NP_002146</u>
Synonyms:	HSP70B'
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_002146</u>
Locus ID:	3310
Cytogenetics:	1q23.3
Protein Pathways:	Antigen processing and presentation, Endocytosis, MAPK signaling pathway, Spliceosome



[View online »](#)