

Product datasheet for **TB407915**

GCLC CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GCLC (NM_001498), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207915
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GCLC Mouse Monoclonal Antibody [Clone ID: OT11A3] (TA507318)
ACCN:	NM_001498 , NP_001489
Synonyms:	GCL; GCS; GLCL; GLCLC
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001489
Locus ID:	2729
Cytogenetics:	6p12.1
Protein Families:	Druggable Genome
Protein Pathways:	Glutathione metabolism, Metabolic pathways



[View online »](#)