

Product datasheet for **TB408134**

Mannose Phosphate Isomerase (MPI) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MPI (NM_002435) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208134
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OT14C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Mannose Phosphate Isomerase (MPI) Mouse Monoclonal Antibody [Clone ID: OT11C7] (TA504749)
ACCN:	NM_002435 , NP_002426
Synonyms:	CDG1B; PMI; PMI1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002426
Locus ID:	4351
Cytogenetics:	15q24.1
Protein Families:	ES Cell Differentiation/IPS
Protein Pathways:	Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Metabolic pathways



[View online »](#)