

Product datasheet for **TB408272**

PSMA2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PSMA2 (NM_002787) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208272
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	PSMA2 Mouse Monoclonal Antibody [Clone ID: OTI3E2] (TA505475)
ACCN:	NM_002787 , NP_002778
Synonyms:	HC3; MU; PMSA2; PSC2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002778
Locus ID:	5683
Cytogenetics:	7p14.1
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Proteasome



[View online »](#)