

Product datasheet for **TB408390**

HDAC8 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HDAC8 (NM_018486), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208390
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HDAC8 Mouse Monoclonal Antibody [Clone ID: OT11A4] (TA809651)
ACCN:	NM_018486 , NP_060956
Synonyms:	CDA07; CDLS5; HD8; HDACL1; KDAC8; MRXS6; RPD3; WTS
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060956
Locus ID:	55869
Cytogenetics:	Xq13.1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)