

## Product datasheet for **TB408463**

### WAVE 1 (WASF1) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of WASF1 (NM_001024934), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining   |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC208463  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| ACCN:                                 | <u><a href="#">NM_001024934</a></u> , <u><a href="#">NP_001020105</a></u>   |
| Synonyms:                             | NEDALVS; SCAR1; WAVE; WAVE1   |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <u><a href="#">NP_001020105</a></u>   |
| Locus ID:                             | 8936  |
| Cytogenetics:                         | 6q21  |
| Protein Families:                     | Druggable Genome  |
| Protein Pathways:                     | Adherens junction, Fc gamma R-mediated phagocytosis, Regulation of actin cytoskeleton   |



[View online »](#)