

Product datasheet for **TB408595**

TRP1 (TYRP1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TYRP1 (NM_000550) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208595
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_000550 , NP_000541
Synonyms:	b-PROTEIN; CAS2; CATB; GP75; OCA3; TRP; TRP1; TYRP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000541
Locus ID:	7306
Cytogenetics:	9p23
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Melanogenesis, Metabolic pathways, Tyrosine metabolism



[View online »](#)