

Product datasheet for **TB408720**

GABPB2 (GABPB1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GABPB1 (NM_005254), transcript variant beta-1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208720
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GABPB2 (GABPB1) Mouse Monoclonal Antibody [Clone ID: OTI7A7] (TA811268)
ACCN:	<u>NM_005254</u> , <u>NP_005245</u>
Synonyms:	BABPB2; E4TF1; E4TF1-47; E4TF1-53; E4TF1B; GABPB; GABPB-1; GABPB2; NRF2B1; NRF2B2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_005245</u>
Locus ID:	2553
Cytogenetics:	15q21.2
Protein Families:	Transcription Factors



[View online »](#)