

Product datasheet for **TB408766**

AMPK beta 2 (PRKAB2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PRKAB2 (NM_005399) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208766
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	AMPK beta 2 (PRKAB2) Mouse Monoclonal Antibody [Clone ID: OTI4H4] (TA808333)
ACCN:	NM_005399 , NP_005390
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005390
Locus ID:	5565
Cytogenetics:	1q21.1
Protein Families:	Druggable Genome
Protein Pathways:	Adipocytokine signaling pathway, Hypertrophic cardiomyopathy (HCM), Insulin signaling pathway



[View online »](#)