

## Product datasheet for **TB408793**

### **C15orf38 (ARPIN) CytoSection**

#### **Product data:**

|  |   |
|--|---|
| <b>Product Type:</b>                         | CytoSections  |
| <b>Description:</b>                          | Transient overexpression of ARPIN (NM_182616) in HEK293T cells paraffin embedded controls for ICC/IHC staining  |
| <b>Species:</b>                              | Human   |
| <b>Expression Host:</b>                      | HEK293T   |
| <b>Expression cDNA Clone or AA Sequence:</b> | TrueORF Clone RC208793  |
| <b>Tag:</b>                                  | C-MYC/DDK   |
| <b>Detection Antibodies:</b>                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)   |
| <b>ACCN:</b>                                 | <u><a href="#">NM_182616</a></u> , <u><a href="#">NP_872422</a></u>   |
| <b>Synonyms:</b>                             | C15orf38  |
| <b>Storage:</b>                              | Room Temperature, or 2-8°C for long term storage  |
| <b>Stability:</b>                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| <b>Preparation:</b>                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| <b>Note:</b>                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| <b>RefSeq:</b>                               | <u><a href="#">NP_872422</a></u>  |
| <b>Locus ID:</b>                             | 348110  |
| <b>Cytogenetics:</b>                         | 15q26.1   |



[View online »](#)