

## Product datasheet for **TB408854**

### H4C3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HIST1H4C (NM_003542) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208854
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_003542</a> , <a href="#">NP_003533</a>
Synonyms:	dj221C16.1; H4-16; H4/g; H4C1; H4C2; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; H4FG; HIST1H4C
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003533</a>
Locus ID:	8364
Cytogenetics:	6p22.2
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)