

## Product datasheet for **TB408974**

### USP38 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of USP38 (NM_032557) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208974
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	USP38 Mouse Monoclonal Antibody [Clone ID: OTI1D11] (TA500312)
ACCN:	<a href="#">NM_032557</a> , <a href="#">NP_115946</a>
Synonyms:	HP43.8KD
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_115946</a>
Locus ID:	84640
Cytogenetics:	4q31.21
Protein Families:	Druggable Genome, Protease



[View online »](#)