

Product datasheet for **TB409446**

LYPLA3 (PLA2G15) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PLA2G15 (NM_012320) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209446
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_012320 , NP_036452
Synonyms:	ACS; GXVPLA2; LLPL; LPLA2; LYPLA3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_036452
Locus ID:	23659
Cytogenetics:	16q22.1
Protein Families:	Secreted Protein
Protein Pathways:	Glycerophospholipid metabolism, Lysosome



[View online »](#)