

## Product datasheet for **TB409510**

### MASA (ENOPH1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ENOPH1 (NM_021204) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209510
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	MASA (ENOPH1) Mouse Monoclonal Antibody [Clone ID: OTI6B6] (TA810832)
ACCN:	<u><a href="#">NM_021204</a></u> , <u><a href="#">NP_067027</a></u>
Synonyms:	E1; MASA; MST145; mtnC
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_067027</a></u>
Locus ID:	58478
Cytogenetics:	4q21.22
Protein Pathways:	Cysteine and methionine metabolism



[View online »](#)