

Product datasheet for **TB409543**

GTF2H2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GTF2H2 (NM_001515) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209543
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_001515</u> , <u>NP_001506</u>
Synonyms:	BTF2; BTF2 p44; BTF2P44; p44; T-BTF2P44; TFIIH
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001506</u>
Locus ID:	2966
Cytogenetics:	5q13.2
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transcription Factors
Protein Pathways:	Basal transcription factors, Nucleotide excision repair



[View online »](#)