

## Product datasheet for **TB409568**

### CD73 (NT5E) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NT5E (NM_002526), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209568
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD73 (NT5E) Mouse Monoclonal Antibody [Clone ID: OT11G2] (TA809085)
ACCN:	<a href="#">NM_002526</a> , <a href="#">NP_002517</a>
Synonyms:	CALJA; CD73; E5NT; eN; eNT; NT; NT5; NTE
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_002517</a>
Locus ID:	4907
Cytogenetics:	6q14.3
Protein Families:	ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism, Purine metabolism, Pyrimidine metabolism



[View online »](#)