

Product datasheet for **TB409595**

HA1 (ARHGAP45) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ARHGAP45 (NM_012292) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209595
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_012292</u> , <u>NP_036424</u>
Synonyms:	HA-1; HLA-HA1; HMHA1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_036424</u>
Locus ID:	23526
Cytogenetics:	19p13.3



[View online »](#)